


HAEMOLANCE

Instructions for Use:

1 Remove Cap


Twist back and forth until loose to easily remove the protective cap.



2 Position Lancet

Place lancet against side of fingertip and press button firmly.

After use, dispose Haemolance Plus in a sharps container.



Product Name	Product Description	Product Number
Haemolance Plus Micro Flow 150 ct.	28 gauge (Light Blue)	990200
Haemolance Plus Micro Flow 50 ct.	28 gauge (Light Blue)	990250
Haemolance Plus Low Flow 150 ct.	25 gauge (Blue)	990900
Haemolance Plus Low Flow 50 ct.	25 gauge (Blue)	990950
Haemolance Plus Normal Flow 150 ct.	21 gauge (Green)	990700
Haemolance Plus Normal Flow 50 ct.	21 gauge (Green)	990750
Haemolance Plus High Flow 150 ct.	18 gauge (Yellow)	990800
Haemolance Plus Pediatric Flow 150 ct.	1.2 mm (Pink)	990400
Haemolance Plus Max Flow 150 ct.	1.6 mm (Purple)	990600

Distributed by:



Minneapolis, MN 55439 USA

800-818-8877

Fax: 800-307-3334

www.hypoguard.com

Made in Sweden

P/N 6070-02

Rev 9/04